

## DETECTION OF ALLERGEN-SPECIFIC IgE

## FIELD OF THE INVENTION

This invention relates to point-of-care test (POCT) devices and methods to screen  
5 animals showing signs of allergy or atopic disease. More particularly, the immunological  
devices and methods of the invention use a defined mixture of clinically relevant  
allergens to detect allergen-specific IgE in samples obtained from such animals.

## CROSS REFERENCE TO RELATED APPLICATIONS

10 This application claims priority to U.S. Provisional Patent Application Serial  
No. 60/325,812, filed on September 28, 2001 and U.S. Provisional Patent Application  
Serial No. 60/259,450, filed on January 3, 2001, both entitled "IN-CLINIC TEST FOR  
THE DETECTION OF ALLERGEN-SPECIFIC IgE" and are incorporated herein by  
reference in their entirety.

15

## BACKGROUND OF THE INVENTION

The cross-linking of mast cell-bound allergen-specific IgE to allergens induces  
Type I allergic diseases such as atopic dermatitis and atopic asthma. Diseases related to  
allergy and atopy affect a significant percentage of the population, including up to 20%  
20 of humans, and are increasing every year.

Allergy-related illness is not restricted to humans, and can be a major problem in  
veterinary practice as well. A number of other diseases may present with symptoms  
similar to those of allergy/atopy, thus a definitive diagnosis may be difficult. A definitive  
diagnosis relies upon a thorough history, careful review of clinical signs, and *in vitro* and  
25 *in vivo* diagnostic tests.

Immunoassay procedures have provided sensitive diagnostic tools for the *in vitro* detection of a variety of allergens of clinical significance. These procedures use allergens and IgE specific detection reagents to measure allergen-specific antibodies in the serum of the patient. These processes characteristically require an incubation period with both the allergens and other reagents and, as a result, may be time consuming. It is not unusual for an assay procedure to require up to 24 hours to complete. Furthermore, the need to adhere to time incubation steps and multiple washings with measured reagents has largely limited these procedures to large hospital and reference clinical laboratories where highly trained personnel and sophisticated equipment are available to perform the assays. Additionally, a certain percentage of atopic patients are allergic to allergens not present in the standard allergen-panels used by reference labs. In these cases, detailed testing fails to reveal the aggravating allergen(s) resulting in these patients incurring unnecessary medical expenses.

Presently there is an unmet need for a simple and rapid procedure for conducting immunoassays for allergen specific IgE which can be used in the physician's or veterinarian's office. The present invention satisfies this need and provides related advantages as well.

#### SUMMARY OF THE INVENTION

The present invention provides rapid and simple tests to aid in the diagnosis of allergy in animals, including, but not limited to, humans, dogs, cats and horses. More particularly, the devices and methods of the invention screen for allergen-specific IgE in a sample using immunological reactions. Such samples can be obtained from animal bodily fluids that contain IgE, such as, for example, urine, mucosal secretions, milk, whole blood or components thereof including plasma, serum and the like.

The devices of the present invention include a first member containing a porous member, such as a membrane or filter, having a sample receiving area. The porous member can be spotted with at least one mixture of allergens (in one spot or multiple spots) and at least one control spot. The control spot or spots can be a positive control, such as purified IgE, and/or a procedural control, such as biotinylated bovine serum albumin (BSA), for example. The mixture of allergens can be spotted either within the sample receiving area or in a separate area of the porous member.

The device can optionally contain an absorbent member in capillary communication with the porous first member to induce the flow of a liquid through the porous member without the use of external means. The absorbent member is selected to have a capillary pore size to induce flow of liquids through the porous member when the hydrostatic pressure of the sample and subsequent reagents are insufficient to induce flow through the porous member.

The device also includes a support member to support or house the first porous member. The absorbent member can also provide the means to support the first member.

Methods for detecting the presence of allergen-specific IgE in a sample use simple devices of the present invention, small volumes, and do not require lengthy incubation steps. Such methods are accomplished by:

- (a) adding the sample to the sample receiving area of a device of the present invention, whereby at least one allergen in the mixture of allergens spotted on the porous member of the device binds to allergen-specific IgE if present in the sample;
- (b) contacting the mixture of allergens with a reagent that selectively binds the allergen-specific IgE; and
- (c) detecting bound allergen-specific IgE if present in the sample.

Useful mixtures of allergens for spotting onto the porous member can contain, but are not limited to, one or more types of allergens selected from the following: cat allergens, dog allergens, flea saliva proteins, house dust mite allergens, storage dust mite allergens, Japanese Cedar pollen, ragweed allergens, ryegrass allergens, meadow fescue, orchard grass, Bermuda grass, sheep sorrel, Kochia, Russian thistle, yellow dock, timothy grass, redtop, birch and blue grass allergens.

Useful reagents that selectively bind the allergen-specific IgE in a sample include, for example, anti-IgE antibodies or IgE receptors or mimetopes of such molecules. Such reagents can be labeled with a detectable marker or a detectable marker can be separately added that binds or otherwise associates with the reagent such that detection is possible.

The present invention accordingly provides low cost point-of-care (i.e., in clinic or field) tests that can be performed simply and quickly in the physician's or veterinarian's office. The rapid results available with the methods and devices of the present invention permit the physician or veterinarian to quickly assess the need for further diagnostic work

and/or treatment plans. Further diagnostic work may include submission of a sample to a central diagnostic laboratory for testing to determine the exact allergens for which the animal is positive, and for which immunotherapy would be appropriate.

Thus, the present invention further provides methods for prescribing immunotherapy treatment of an animal having an allergic disease. Briefly, the methods are accomplished by:

- (a) obtaining a sample from the animal;
- (b) analyzing the sample according to the methods of the present invention to detect the presence of allergen-specific IgE in the sample, wherein the presence of allergen-specific IgE identifies the sample as a positive sample;
- (c) analyzing the positive sample to identify one or more allergens responsible for an allergic disease; and
- (d) selecting an appropriate immunotherapy to treat the allergic disease.

The present invention further provides novel reagents, including anti-IgE antibodies, and kits containing devices and ancillary reagents for conducting the immunoassays.

## DETAILED DESCRIPTION OF THE INVENTION

The present invention is generally directed to POCT devices and immunoassay methods for the rapid detection of allergen-specific IgE. The devices and methods of the present invention are useful as preliminary aids in the diagnosis of allergy in animals.

The immunoassay methods of the present invention use POCT devices that include a first porous member and a support member. Particularly suitable devices include the following elements: (1) a porous solid surface to which a defined mixture of allergens can be spotted in one or more spots; (2) a sample receiving area for adding the sample to be assayed; (3) a means for drawing liquids through the porous surface; and (4) an area for displaying the results. In certain embodiments, one element may perform the functions of one or more of the other elements. For example, the sample receiving area may be the same area or a different area from where the defined mixture of allergens have been spotted.

The porous first member can be made of any suitable material to which the mixture of allergens can be bound or affixed, while sufficiently porous to allow the flow of liquids to pass through or by the affixed mixture of allergens. Suitable materials include membranes, filters and the like. Such membranes and filters can be made from a variety of materials, which include but are not limited to, for example, plastic, glass, gel, celluloid, paper, fabric, polystyrene, nylon, nitrocellulose, agarose, cotton, and PVDF (poly-vinylidene-fluoride). The mixture of allergens can be spotted (i.e. affixed) on the porous member by any means known to those skilled in the art or as described in the Examples below.

One or more controls can also be spotted on the porous member to ensure the proper operation of the immunoassay. Preferably, the control spot is an area on the device, which comes into contact with the assay reagents and which contains an agent which serves to indicate binding of IgE-detection reagent has occurred. Agents useful as positive controls include, for example, purified polyclonal or monoclonal IgE, IgE receptors or active fragments thereof, anti-IgE immunoglobulin, anti-IgG immunoglobulin, protein A, mimetopes of any of the above molecules or any other molecule resulting in binding of the IgE-detecting reagent to the porous member. As used herein, the term "mimeto" refers to any compound able to mimic the ability of anti-IgE antibody to bind to IgE. An example of a mimeto can be a peptide that has been modified to decrease its susceptibility to degradation but which still retains the ability to bind the IgE-detecting reagent. Other examples of mimetopes include, but are not limited to, carbohydrate-based compounds, lipid-based compounds, nucleic acid-based compounds, natural organic compounds, synthetically derived organic compounds, anti-idiotypic antibodies and/or catalytic antibodies, or fragments thereof. A mimeto can be obtained by, for example, screening libraries of synthetic compounds for compounds capable of binding to the IgE-detecting reagent. Procedural controls known to those skilled in the art can also be spotted on the porous member, including, for example, biotinylated BSA. The porous first member further contains

a sample receiving area. As noted previously, the sample receiving area can be the same area where the mixture of allergens are spotted or can be in a different location on the porous member depending on the desired format. A pre-filter can optionally be used for samples that are preferably filtered prior to introducing the filtered sample to the sample receiving area, including, for example, whole blood to obtain serum. Such samples can be filtered separately from the device of the present invention in a pre-filter device or the pre-filter device can be seated above the sample receiving area and, if desired, the pre-filter can be removed and disposed after the sample has been pre-filtered.

The devices of the present invention optionally contain an absorbent member having capillary pathways connected to the porous member to permit a liquid to traverse the absorbent member. Accordingly, the absorbent member is in capillary communication with the porous member and has capillary pores sized to induce the flow of liquid through the first porous member without the use of external means. The absorbent member is particularly useful when the hydrostatic pressure of the sample and subsequent reagents used in the assay are insufficient to induce flow through the first member. In addition, the absorbent member can also provide support for the first member.

Those skilled in the art can readily identify suitable materials useful as absorbent members and other support members. Thus, the devices of the present invention can be of any standard format known to those skilled in the art, including, for example, a lateral flow device (as described in U.S. Patent 5,073,484 or U.S. Patent 5,654,162, both incorporated herein by reference), a flow-through device (as described in U.S. Patent 4,727,019, incorporated herein by reference), or any modifications of such devices readily known or determined by those skilled in the art.

The methods of the present invention for detecting the presence of allergen-specific IgE in a sample obtained from an animal are generally accomplished by first adding a sample to the sample receiving area of a device as described above. After the sample has come in

contact with the test spot where the mixture of allergens has been spotted on the porous member, a reagent that selectively binds allergen-specific IgE is added and allowed to come in contact with the test spot. If present in the sample, the bound allergen-specific IgE is then detected.

5 As used herein, the term "animal" refers to any mammal that is capable of having an allergic disease caused by the presence of deleterious amounts of allergen-specific IgE. Such animals preferably include, but are not limited to, humans, dogs, cats and horses. As used herein, a cat refers to any member of the cat family (i.e., Felidae), including domestic cats, wild cats and zoo cats. Examples of cats include, but are not limited to, domestic cats, lions, 10 tigers, leopards, panthers, cougars, bobcats, lynx, jaguars, cheetahs, and servals. A preferred cat to test is a domestic cat. As used herein, a dog refers to any member of the family Canidae, including, but not limited to, domestic dogs, wild dogs, foxes, wolves, jackals, and coyotes and other members of the family Canidae. As used herein, a horse refers to an equid. An equid is a hoofed mammal and includes, but is not limited to, domestic horses and wild 15 horses, such as, horses, asses, donkeys, and zebras. Preferred horses to test include domestic horses, including race horses.

Samples to be assayed can be obtained from bodily fluids of such animals. Any bodily fluid that contains allergen-specific IgE can be used, including, for example, urine, mucosal secretions, milk, whole blood and components thereof, such as serum and plasma. 20 Further examples of bodily fluids may be found in U.S. Patent 6, 309, 832 which is hereby incorporated by reference. A bodily fluid can be used undiluted as the sample or can be diluted as desired in a suitable diluent before applying the sample to the sample receiving area. Those skilled in the art can readily determine the appropriate sample to be used in a particular assay.

25 The allergens that are selected in the defined mixture of allergens that is deposited on the porous member are selected such that a minimum number of spots, preferably one spot,

will identify a great majority of patients with IgE-mediated allergic disease. Preferably, the defined mixture of allergens will be capable of identifying at least 50%, 55%, 60%, 65%, 70%, 75%, 80%, 85%, 90%, 95%, 99% and 100% of patients with IgE-mediated allergic disease. Those skilled in the art can readily determine a defined mixture of allergens capable

5 of identifying a desired percent of the relevant population. For example, if the desired target is 90% of all human patients with IgE-mediated allergic disease, those skilled in the art would select a mixture of allergens that are collectively known to cause 90% of all IgE-mediated allergies in humans. Ideally, the mixture of allergens is chosen such that the minimum number of allergens allows the identification of the maximum number of patients.

10 Various embodiments of the present invention may have different mixtures of allergens based on the patient population to be tested.

The allergens may be in the form of extracts, purified allergenic proteins from extracts or from raw materials such as mites or pollens, recombinant allergens, or combinations of the above. Suitable allergens for use in the defined mixture include, but are

15 not limited to, cat allergens (U.S. Patent 5,328,991; 6,019,972; 6,048,962; and 6,120,769); dog allergens (U.S. Patent 5,843,672); flea saliva proteins (U.S. Patent 5,646,115; 5,795,862; and 5,840,695); house dust mite allergens (U.S. Patent 5,433,948; 5,552,142; 5,770,202; 5,773,002; 5,820,862; 5,986,526; 5,972,352; 6,060,057; 6,071,522; 6,077,518; 6,086,897; 6,132,734; and 6,077,517); Japanese Cedar pollen (U.S. Patent 6,090,386); Johnson grass

20 allergens (U.S. Patent 5,480,972; 5,691,167; and 5,736,149); ragweed allergens (U.S. Patent 5,698,204); ryegrass allergens (U.S. Patent 5,710,126; 5,840,316; and 5,965,455); meadow fescue, orchard grass, Bermuda grass, sheep sorrel, Kochia, Russian thistle, yellow dock, timothy grass, redtop, and blue grass allergens. Additional allergens may be found in U.S. Patent 6, 309, 832. The above U.S. patents are hereby incorporated by reference.

25 In the methods of the present invention, after the sample comes in contact with the spotted mixture of allergens in order to allow any allergen-specific IgE to bind to one or



more of the allergens in the mixture, a reagent that selectively binds allergen-specific IgE is then added and allowed to bind to any bound allergen-specific IgE. The term "second binding reagent" is used interchangeably with "reagent that selectively binds allergen-specific IgE." The detection of allergen-specific IgE bound to the second binding reagent is  
 5 indicative of the presence of allergen-specific IgE in the sample.

As used herein, a reagent that "selectively binds allergen-specific IgE" means that it preferentially binds IgE as compared to IgG or another class of immunoglobulin such that the second binding reagent does not bind to IgG or another class of immunoglobulin in an amount that substantially affects the results of the immunoassay. For example, a suitable  
 10 second binding agent binds IgE to a significantly greater extent (has a substantially higher affinity for IgE) than IgG. In a preferred embodiment, the second binding agent binds IgE with a specificity similar to that exhibited by the human IgE receptor alpha chain protein. In another embodiment, the second binding agent fails to bind IgG under conditions in which known anti-IgE antibodies bind IgE. Binding of immunoglobulin by binding agents of the  
 15 present invention can be measured using procedures known to those skilled in the art, for example, the ELISA procedures described herein. In a particularly preferred embodiment, the second binding agent binds IgG with an absorbance value, as determined by ELISA, of less than 300, 250, 225, 200, 175, 150, 125 or 100. Guidance for selecting suitable binding agents is provided in the Examples below.

20 The second binding reagent may be a polyclonal or monoclonal antibody that selectively binds IgE, the Fcε receptor 1 (FcεR1) (as described in U.S. Patent 5,945,294, which is incorporated herein by reference), IgE-binding allergens or any molecule which selectively binds to IgE. One example of such an antibody is an anti-IgE antibody that only recognizes antibodies which bind to mast cells and basophils. Preferred antibodies of the  
 25 present invention have the characteristics of a mAb produced using human FcεRIα-affinity purified IgE as described herein.

In a preferred embodiment, the second binding reagent is a monoclonal antibody that selectively binds to an IgE receptor-binding region or domain. Examples of such antibodies include, but are not limited to, for example, monoclonal antibodies H-219, H-195, H-202, H-203 and H-207. Species-specific antibodies and receptors are particularly useful for

5 detecting allergen-specific IgE in samples obtained from that species. For example, canine anti-IgE antibodies are useful as the second binding reagent in assays for detecting allergen-specific IgE in sample obtained from a dog

In a preferred embodiment, the second binding reagent is monoclonal antibody that selectively binds to canine IgE. Examples of such antibodies include, but are not limited to,

10 for example, monoclonal antibodies H-219, H-220, H-189, H-190, H-195, H-197, H-200, H-202, H-203, H-206 and H-207. In another preferred embodiment, the second binding reagent is monoclonal antibody that selectively binds to feline IgE. Examples of such antibodies include, but are not limited to, for example, monoclonal antibodies H-99, H-100, H-101, H-102, H-106. In another preferred embodiment, the second binding reagent is monoclonal

15 antibody that binds to both canine and feline IgE. Examples of such antibodies include, but are not limited to, for example, monoclonal antibodies H-201, H-192, H-198, H-199, H-208. In another preferred embodiment, the second binding reagent is monoclonal antibody that selectively binds to equine IgE. Examples of such antibodies include, but are not limited to, for example, monoclonal antibody H-357.

20 In a further preferred embodiment, the second binding reagent is an IgE receptor. The term "IgE receptor" includes the full length receptor and active fragment thereof. As used herein, the term "active fragment thereof" means fragments of the full length receptor having binding sites that are capable of binding to the allergen-specific IgE captured by the allergen(s) contained in the defined mixture.

25 In one embodiment, the second binding reagent can be directly conjugated to, or labeled with, an indicator agent. In an alternative embodiment, an indicator agent capable of

binding to or associating with the second binding reagent can be separately added and detected. Those skilled in the art can readily identify suitable detection means for use in the present methods.

Any indicator or detection reagent such as enzymes, fluorescent molecules,

5 radiolabels, or other detectable marker known to those skilled in the art can be used.

Enzymes are particularly useful as indicator agents. If the label is an enzyme, after washing to remove unbound receptor, a solution of the enzyme substrate is added. If the device has bound to it the allergen specific IgE, then the IgE will have bound to it a portion of labeled receptor. The enzyme will cause the substrate to react and generate, if properly selected, a  
10 new, and easily detectable enzymatic-reaction end product. A particularly useful end product is a chromophore which results in a visual color change. If a label other than an enzyme is used, the procedure can be varied. Fluorescence on the porous member can be measured by known methods if a fluorescing label is used. If a radionuclide label, such as  $^{125}\text{I}$  is used, the porous member can be removed and counted.

15 The addition of the sample, second binding reagent and any desired washing buffers can be preceded by brief incubation periods to permit more extensive binding of the allergen-specific IgE in the sample to increase the sensitivity of the assay. However, in the preferred methods of the present invention, such incubation steps are either unnecessary or are very brief, i.e., generally on the order of 60 seconds or less. The flow of solutions containing the  
20 sample or second binding reagent results in a substantially faster rate of binding than is observed in the absence of flow. Accordingly, the total time required for performing the test is about 5 minutes or less, and preferably less than about 3 minutes.

The present invention also provides methods of prescribing immunotherapy treatment in animals having IgE-mediated allergic disease. Such methods include obtaining a sample  
25 from an animal having symptoms of allergic disease, analyzing the sample according to the immunoassays methods described above to detect the presence of allergen-specific IgE. If

allergen-specific IgE is detected in the sample, the positive sample is then further analyzed to identify one or more allergens responsible for the allergic disease and an appropriate immunotherapy is then selected for the animal. Further analysis of a positive sample to identify the specific allergen or allergens responsible for the allergic disease can be

5 accomplished by any means known to those skilled in the art, including, for example, by *in vitro* immunoassays conducted for example in a reference laboratory or by intradermal tests (e.g. conventional skin tests). Once the specific allergen or allergens have been identified, clinicians skilled in the art can determine the appropriate immunotherapy for the particular patient.

10 Kits for performing the methods of the present invention are also provided. In addition to the devices of the present invention, the kits of the present invention include reagents used in the methods of the present invention. Such reagents include wash buffers and detection reagents. Preferred detection reagents of the invention include monoclonal antibodies that specifically bind canine, feline, equine, and human IgE.

15 The following examples are provided for the purposes of illustration, and are not intended to limit the scope of the present invention.

### Example 1

20 This example describes how the allergens for the screens are selected, and how the devices are prepared.

The methods of the present invention use carefully defined mixtures of allergens. The mixtures are selected to identify most allergic individuals with a minimum allergen panel.

25 Data from the assignee's US diagnostic laboratories was analyzed to identify allergens to which many dogs were allergic. The point of care test was designed to detect 70 to 90% of animals that give positive results in the complete panel at the diagnostic laboratories, while keeping the total number of allergens in the point of care test small.

A total of more than two thousand samples was analyzed in several different sets. The samples were tested in the full panel of allergens and approximately 80% tested positive to at

least one allergen in the panel. (An animal that tests negative to the complete panel may be an animal that is allergic to an allergen not in the panel, or may not be allergic at all, but rather have another disease with symptoms that resemble allergy).

Of those 80%, approximately 60% are positive to either flea or dust mites. To add an additional 15% of animals, other important allergens may be added. Candidate additional allergens would be Meadow Fescue, Orchard Grass, red Cedar, yellow dock, timothy grass, redtop, and blue grass. A different set of allergens may be desired in different geographic regions and/or climates, and/or for different animals (the above data is for dogs). The desired set of allergens is again determined by analysis of allergen specific IgE detected in a standard or reference laboratory tests.

Data from other regions was analyzed in a similar manner. Nearly 800 dog sera were analyzed for antibodies to indoor allergens. House dust mite was the most common allergen, followed by storage mites and flea saliva. A point of care device with house dust mite, storage mite, and flea saliva could detect 98% of the animals that are positive in the indoor panel.

More than 500 dog sera were analyzed for United Kingdom general pollens; 72 % of the dogs were positive to cocksfoot. The next most prevalent pollen allergens were alder, English Plantain, Lamb's quarters, and Sheep Sorrel. A point of care device with allergens from house dust mite, storage mite, and cocksfoot or highly related grasses, could detect 91% of the animals that are positive in the complete panel.

Approximately 30 dog sera were analyzed for Central Europe panel; 100 % of the dogs were positive to house dust mite or storage mites. A point of care device with house dust mite and storage mite could detect 100% of the animals that are positive in the complete panel.

Approximately 56 dog sera were analyzed for Northern Europe panel; 65 % of the dogs were positive to house dust mite. The next most prevalent allergens were storage mites, followed by red top and velvet grass. A point of care device with house dust mite and storage mite could detect 98% of the animals that are positive in the complete panel.

Approximately 122 dog sera were analyzed for United Kingdom extended pollens; 58 % of the dogs were positive to Meadow fescue. The next most prevalent allergen was Grand

5      fir. A point of care device with Meadow fescue could detect 58% of the animals that are positive in the complete panel.

Approximately 104 dog sera were analyzed for Mediterranean panel; 59 % of the dogs were positive to house dust mite. The next most prevalent allergens were storage mites and flea saliva. A point of care device with house dust mite, storage mite, and flea saliva could detect 97% of the animals that are positive in the complete panel.

To determine the correct allergens to use in the devices and methods of the invention it is also necessary to assess such allergens in a point of care format, to determine that the sera will still test positive when the allergens are present as a mixture, rather than as an isolated sample.

In one embodiment of the present invention, the point of care device is prepared by spotting the porous member with one mixture of allergens comprising meadow fescue, dust mite, and flea saliva allergens. Each of these allergens may be present as crude extract, purified extract, and/or genetically engineered and recombinantly expressed pure proteins. The allergens extracts or recombinant proteins are spotted onto the porous member at approximately 1,000 to 500,000 PNU/ml; or 10 µg/ml to 10 mg/ml, depending upon the particular allergen. The exact concentration required for a test must be determined experimentally, but one can begin with a concentration similar to that used in an ELISA format with the same allergen or extract. The allergens are spotted onto the porous member with spotting buffer without BSA (TRIS 6.05 g/l; NaCl 8.5 g/l; NaN<sub>3</sub> 20 mg/l; and sucrose 10 g/l).

A spot of canine IgE is added at a different location as positive control. This is spotted at 1 µg/ml to 10 mg/ml, preferably about 5-10 µg/ml in spotting buffer with BSA (TRIS 6.05 g/l; NaCl 8.5 g/l, NaN<sub>3</sub> 20 mg/l, sucrose 10 g/l, and BSA 1g/l).

In another embodiment of the invention, there are two test spots on the device, one with meadow fescue and house dust mite, and the other with flea saliva proteins, in addition to the control spot of canine IgE.

Example 2

This example describes the production of mouse anti-canine IgE monoclonal antibodies (mAb) essentially by the method of Kohler and Milstein (Kohler and Milstein, (1975), Nature 256:495-497), which is herein incorporated by reference.

Canine IgE was prepared by pooling canine serum and fractionating by ammonium sulfate precipitation. The fraction from 30-55% saturation was collected and passed through a Protein G column. The flow through fraction was applied to a rHuman IgE Receptor Affiprep (Biorad) column. Retained IgE was recovered from the column with high pH elution. A polybuffer Exchanger column further purified IgE. Fractions eluting at pH5.7 to 4.95 were pooled.

Two Balb/c mice were immunized in the footpad with 30 µg canine IgE suspended in phosphate buffered saline (PBS) and Freund's complete adjuvant. A boost of 30 µg canine IgE was given in PBS/Freund's incomplete adjuvant in the footpad on day 14. Sera was tested for presence of anti-canine IgE antibodies at day 21. The mouse exhibiting the highest titer against canine IgE by ELISA was boosted with 5 µg canine IgE in PBS, intravenously on day 42. Spleenocytes were harvested three days later, and fused with mouse SP2/0 myeloma cells at mid-log growth phase using PEG.

Cells were cultured in RPMI medium containing 20% fetal bovine serum, 10% thymocyte conditioned media, 2 mM L-glutamine, 1 mM sodium pyruvate, and 60 µM β-mercaptoethanol. Hybrids were selected by adding 100 µM hypoxanthine, 10 µM thymidine and 0.4 µM aminopterin. Wells containing hybridoma colonies were tested for specific mAb production against canine IgE, IgG, and IgM by ELISA. Antigens were coated at 5 ng/well. Wells positive to canine IgE were expanded and dilution cloned by limiting dilution until stable mAb secreting cell lines were identified.

Monoclonal antibody is either supplied in a crude form as tissue culture supernatant, ascites or bioreactor supernatant. Specific mAb can be purified using affinity chromatography or ion exchange chromatography from any of these fluid types.

Twenty antibodies reactive to canine IgE epsilon chain by ELISA were generated. All mAbs recognize canine IgE and not canine IgG or IgM by ELISA. Of these, fifteen recognize the 71 kDa epsilon chain by Western blots except three (H-198, 199, and 206) and two are yet to be tested (H-192 and 204). Three Western blots patterns emerged using these

mAbs to probe canine IgE; one set of mAbs (H-195, 197, and 208) recognize the 71 kDa epsilon chain of canine IgE only; a second set of mAbs (H-189, 190, 191, 193, 203, 207, 219 and 220) recognize the 71 kDa epsilon chain of canine IgE and a band running right above the dye front (at about 6 kDa); and the third set of mAbs (H-200, 201, 202, and 205) recognize the 71 kDa epsilon chain of canine IgE and a broad diffuse band at 55 kDa. This third set of mAbs do not recognize canine IgG gamma chain, suggesting that the extra band is not a contaminant of dog IgG. This interesting binding domain appears to have some relationship with the canine IgE epsilon chain. All but one mAb (H-204) were evaluated, and five of these mAbs (H-191, 192, 198, 199, and 208) recognized feline IgE by ELISA even better than canine IgE, and in some cases reactivity with canine IgE was not even detectable. One mAb (H-201) recognized feline IgE equally well as canine IgE.

Antibodies were grouped into 9 categories based on ELISA reactivity patterns. The highest titered mAb from each group was purified and screened as a canine IgE detection POCT candidate (H-189, 199, 200, 201, 204, 205, 206, 207, and 208). From these mAbs, H-200, 201, and 207 were chosen. Further testing with these three mAbs indicated that H-207 was the preferred mAb to be utilized in the POCT.

The antibodies mentioned above can be utilized in various assay configurations to detect various levels of canine IgE present in complex samples such as dog sera, PBMCs or possibly even blood. They work well in combination with the human FcεRIα to calculate total serum IgE. With appropriate standards, dog IgE levels could easily be measured and compared to various stages of immune stimulation, allergen exposure, etc. Additional applications of anti-canine IgE mAbs include application to simplifying purification of canine IgE and/or canine IgE immune complexes by affinity chromatography methods. Use of mAbs specific to epsilon chain would easily and more readily remove dog IgE from sample sources such as serum as compared to conventional chromatography methods. Other uses include monitoring separation of intact canine IgE, epsilon chain during purification processes. Flow cytometry studies could be performed by using these antibodies to map cell surface IgE and to check intact canine IgE, epsilon chain during purification processes. Flow cytometry studies could be performed by using these antibodies to map cell surface IgE and to check the validity of FcεRIα based results. In addition, the anti-canine epsilon chain mAbs may be used for immunotherapy for allergic dogs.



Mouse anti-feline IgE mAbs, mouse anti-equine IgE mAbs, and mouse anti-human IgE mAbs are made in the same manner as above. The mAbs may also be prepared in animals other than mouse by standard techniques for preparation of monoclonal antibodies for the above cited uses.

5 Example 3

This Example discloses the ability of the anti-canine IgE monoclonal antibodies to bind antigen-specific, canine or feline IgE, as measured by solid phase ELISA, using the following protocol:

The wells of a microtiter plate were coated with either Der F protein (30 Protein  
10 Nitrogen Units (PNU)/well) (Center Labs, Port Washington, N.Y.), flea saliva antigen (FSA) (100 ng/well) (prepared as described in U.S. Patent 5,646,115) or canine IgG (1µg/well)(Jackson ImmunoResearch Labs, West Grove, PA), diluted in CBC buffer (50 mM carbonate, pH 9.6) and the plate stored, covered, overnight at 4°C. (Wells coated with IgG serve to measure the level of non-specific immunoglobulin binding.) The following day,  
15 excess fluid was removed from the wells, the plates blotted dry, and 200 µl of Assay buffer (4% fetal calf serum in PBS containing 0.05% Tween-20) were added to each well. Following a 60 minute incubation at room temperature (RT), the wells were washed four times using Wash buffer (PBS containing 0.05% Tween-20), and 100 µl of canine or feline serum (diluted 1:20 in Assay buffer) were added to each well. Some wells (labeled 'Heated  
20 Serum') received serum which had been heated at 56°C for one hour prior to use. ( This procedure is known to denature the IgE-Fc receptor binding region which allows one to determine which antibodies bind IgE in or near this region) The plates were incubated overnight at 4°C, after which the wells were washed four times using Wash buffer, and 100 µl of supernatant from anti-canine IgE hybridoma cells (prepared as described in Example 2),  
25 diluted in Assay buffer as indicated, were added to each well. The plates were incubated for 2 hours at room temperature after which, the wells were washed four times using Wash buffer and 100 µl of HRP-conjugated, donkey, anti-mouse IgG (Jackson ImmunoResearch Labs, West Grove, PA), diluted 1:2000 in Assay buffer, were added to each well. After incubating at RT for one hour, the wells were washed four times using Wash buffer and 100  
30 µl of TMB Peroxidase Substrate System (KPL Labs, Gaithersburg, MD) were added to each well. The plates were incubated at RT until color developed within the wells (approximately

- ten minutes), at which time the reaction was stopped by adding 100  $\mu$ l of 1.0M phosphoric acid to each well. Background levels of non-specific binding were measured by preparing wells as described above but without adding canine or feline serum. Using an ELISA plate reader, the optical density of each well was measured at 450 nM, the results of which are
- 5 shown in Table 1. The background level of absorbance was subtracted from all of the values shown in Table 1. Values over 100 are considered to be positive.

Table 1.

mAb ID	Dog Serum				Cat Serum			Dog IgG Serum
	MAb Dilution	Serum	Heated Serum	Serum (-)	Serum	Heated Serum	Serum (-)	
H-219	1:250	3249	479	40	5	22	5	455
	1:500	3028	369		4	4		
	1:1000	2819	270		10	10		
H-220	1:250	49	1533	0	7	6	4	316
	1:500	42	1234		8	4		
	1:1000	60	870		7	6		
H189	1:250	39	1838	2	7	6	8	473
	1:500	37	1714		3	1		
	1:1000	42	1465		2	3		
H-190	1:250	2783	2970	0	65	21	12	1192
	1:500	2227	2798		34	10		
	1:1000	1578	2426		24	7		
H-191	1:250	35	57	0	3630	2636	0	22
	1:500	46	55		3517	2351		
	1:1000	44	60		3245	1714		
H-192	1:250	55	62	8	3418	2926	35	64
	1:500	44	52		3071	2273		
	1:1000	48	71		2595	1808		
H-193	1:250	72	1572	0	13	7	9	368
	1:500	61	1306		11	14		
	1:1000	71	972		8	13		
H-195	1:250	3752	2027	0	26	30	22	1260
	1:500	3604	1477		18	21		
	1:1000	3424	974		21	27		
H-197	1:250	1953	2882	0	28	12	6	902
	1:500	1294	2415		17	11		
	1:1000	739	1779		14	16		
H-198	1:250	120	60	0	3818	2923	1	14
	1:500	81	44		3803	2841		
	1:1000	79	45		3714	2779		
H-199	1:250	60	50	1	3108	2502	0	85
	1:500	44	47		2977	2242		
	1:1000	43	47		2736	1959		
H-200	1:250	3781	2420	4	29	33	28	1355
	1:500	3744	1897		15	17		
	1:1000	3532	1436		11	11		

H-201	1:250	3606	2538	1	3695	3061	26	1782
	1:500	3473	2195		3557	2480		
	1:1000	3223	1904		3446	2705		
H-202	1:250	3697	1892	0	20	17	19	1170
	1:500	3826	1744		15	11		
	1:1000	3620	1343		22	9		
H-203	1:250	3803	2425	0	69	69	68	1986
	1:500	3603	2024		39	42		
	1:1000	3424	1478		32	30		
H-205	1:250	67	3483	0	12	6	9	472
	1:500	66	1245		6	3		
	1:1000	68	752		8	10		
H-206	1:250	1972	1446	0	3	18	8	472
	1:500	1615	1284		4	3		
	1:1000	1156	972		9	11		
H-207	1:250	3393	720	4	12	11	12	663
	1:500	3085	546		8	6		
	1:1000	2617	353		11	10		
H-208	1:250	105	51	1	3053	2939	5	22
	1:500	101	67		2672	2797		
	1:1000	122	59		2184	2488		

#### Example 4

The ability of the anti-canine IgE monoclonal antibodies to bind canine or feline IgE bound by immobilized human Fc-epsilon receptor I-alpha chain, was tested using the following capture-assay protocol.

The wells of a microtiter plate were coated with the extracellular domain of human Fc-epsilon receptor I-alpha chain protein(huFcεRIα) (prepared as described in U.S. patent 5,945,294 ) (42.5 ng/well) in CBC buffer, and the plate incubated, covered, overnight at 4°C. Additionally, several wells were also coated with heat-treated canine IgG (Jackson ImmunoResearch Labs, West Grove, PA) (1 µg/well) to measure the background level of non-specific immunoglobulin binding. It is believed that heat treatment specifically destroys the activity of contaminating IgE.

The following day, excess fluid was removed from the wells, 200 µl of Assay buffer added to each well and the plate incubated at RT for one hour. The wells were washed four times using Wash buffer, and purified canine or feline IgE (prepared as described in Example 2) was added to each well (50 ng/well in 100 µl Assay buffer). Following a two hour incubation at RT, the wells were washed four times using Wash buffer and 100 µl of anti-canine IgE hybridoma cell supernatant (prepared as described in Example 2), diluted as indicated in Assay buffer, were added to each well. The plate was incubated at RT for one

hour, the wells washed four times using Wash buffer, and 100 µl of HRP-conjugated, donkey, anti-mouse IgG diluted 1:2000 in Assay buffer, were added to each well. The plate was incubated for one hour at RT after which, the wells were washed four times using Wash buffer, and 100 µl of TMB Peroxidase Substrate System (KPL Labs, Gaithersburg, MD) were added to each well. The plate was incubated at RT until color developed within the wells (approximately 15 minutes) at which time the reaction was stopped by adding 100 µl of 1.0M phosphoric acid to each well. The optical density of each well was measured at 450 nM using an ELISA plate reader, the results of which are shown in Table 2. Background absorbance was measured by preparing some wells as described above but without any mAb's and background levels were subtracted from all values listed in Table 1. Values over 100 are considered positive.

Table 2.

	H-219			H-220			H-189		
mAb Dilution	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG
1:100	2897	11	77	26	5	1	30	4	2
1:300	2960	6	38	11	5	0	22	0	0
1:1000	2679	5	17	8	0	0	17	4	0
	H-190			H-191			H-192		
mAb Dilution	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG
1:100	3985	152	108	15	576	0	2726	4200	0
1:300	3062	53	67	10	353	0	1199	4200	10
1:1000	2444	17	22	3	137	0	759	3713	6
	H-193			H-195			H-197		
mAb Dilution	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG
1:100	24	5	4	4000	16	127	2570	78	67
1:300	16	5	0	2851	11	75	3132	32	35
1:1000	9	2	0	2069	5	25	2789	14	14
	H-198			H-199			H-200		
mAb Dilution	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG
1:100	2871	4200	3	1240	3181	0	3371	35	137
1:300	1510	4200	71	1091	3344	1	2434	19	120
1:1000	1187	3865	0	993	3068	0	2857	6	41
	H-201			H-202			H-203		
mAb Dilution	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG
1:100	2750	3423	82	3523	68	107	2926	15	177
1:300	3197	3732	88	2938	5	85	3509	22	197
1:1000	3387	3627	71	2000	3	33	3704	13	168
	H-205			H-206			H-207		
mAb Dilution	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG	Dog IgE	Cat IgE	Dog IgG

1:100	27	0	0	2153	5	19	2690	6	35
1:300	6	3	2	2011	6	14	1925	10	16
1:1000	5	0	0	1553	3	4	1414	0	4
H-208									
mAb Dilution	Dog IgE	Cat IgE	Dog IgG						
1:100	131	1880	0						
1:300	219	2099	0						
1:1000	151	1617	64						

### Example 5

This example discloses the ability of anti-feline IgE monoclonal antibodies to bind antigen-specific canine or feline IgE. Hybridoma cell lines producing anti-feline IgE were generated using the procedure described in Example 2, with the exception the mice were immunized with feline IgE instead of canine IgE. The ELISA's were performed using the solid phase ELISA protocol outlined in Example 3 and substituting the anti-feline IgE hybridoma supernatants for anti-canine IgE hybridoma supernatants. Additionally, to measure the background level of non-specific immunoglobulin binding, some wells were coated with feline IgG instead of canine IgG. The results of these ELISA's are shown in Table 3. Background levels of absorbance have been subtracted from all values shown in the table. Values greater than 100 are considered positive.

Table 3.

		Dog Serum				Cat Serum				
		DerF		Flea Saliva		DerF		Flea Saliva		Pure Cat IgG
		Serum	Heated Serum	Serum	Heated Serum	Serum	Heated Serum	Serum	Heated Serum	Serum
H-100	1:250	27	43	40	38	2473	2794	2925	2371	41
	1:500	30	46	39	41	1932	2269	2605	2083	72
	1:1000	39	50	39	42	1287	1582	2001	1479	17
H-101	1:250	56	67	46	45	3577	3464	3741	3287	59
	1:500	41	63	38	47	3144	3042	3392	2966	84
	1:1000	74	61	39	45	2474	2465	2921	2416	21
H-102	1:250	60	54	32	43	3586	3444	3628	3239	33
	1:500	34	46	27	36	3188	3052	3383	2984	53
	1:1000	30	47	29	39	2607	2501	2935	2424	28
H-103	1:250	83	64	43	44	30	35	6	25	0
	1:500	79	72	53	42	3	68	0	5	0
	1:1000	48	54	41	48	7	14	1	7	0
H-106	1:250	53	59	43	43	1268	1443	1580	1002	0
	1:500	42	74	38	50	861	1008	1257	713	5

	1:1000	58	52	45	63	450	537	705	433	0
H-99	1:250	110	205	215	261	3567	4097	3632	3522	4200
	1:500	90	113	130	152	3363	3962	3277	3315	4200
	1:1000	74	107	111	125	3016	3814	3305	2978	4089
H-104	1:250	4200	4200	4200	4200	4200	4200	3895	3929	4200
	1:500	4200	4200	4200	4200	4121	4200	3755	3883	4200
	1:1000	4200	4200	4067	4093	3932	3919	3717	3746	3892
H-107	1:250	49	51	30	42	1418	2226	217	587	3975
	1:500	45	62	15	47	1225	1960	172	504	3808
	1:1000	37	76	40	44	879	1509	143	420	3140

This example demonstrates that all of these mAbs recognize feline IgE. Five out of the eight mAbs are specific for feline IgE whereas two mAb recognizes both feline IgE and IgG and one mAb binds feline IGE, IgG and canine IgE.

### Example 6

This example discloses the ability of the anti-feline IgE monoclonal antibodies to bind huFcεRIα-immobilized, canine or feline IgE. This example uses the capture-assay protocol described in Example 4 with the exception that anti-feline IgE hybridoma cell supernatants were used in place of anti-canine IgE hybridoma cell supernatants. The results of this assay are shown below in Table 4.

Table 4.

	H-99		H-100		H-101	
mAb Dilution	Canine IgE	Feline IgE	Canine IgE	Feline IgE	Canine IgE	Feline IgE
1:00	0	206	24	2496	8	2443
1:300	0	186	11	1591	0	1381
1:1000	0	60	11	846	0	1111
	H-102		H-103		H-104	
mAb Dilution	Canine IgE	Feline IgE	Canine IgE	Feline IgE	Canine IgE	Feline IgE
1:00	3	2325	0	0	333	476
1:300	0	1002	0	0	229	352
1:1000	6	1087	0	0	287	169
	H-106		H-107			
mAb	Canine	Feline	Canine	Feline		

Dilution	IgE	IgE	IgE	IgE		
1:00	9	1238	0	37		
1:300	0	491	0	33		
1:1000	5	245	2	37		

### Example 7

This example discloses a method for the purification of recombinantly produced equine Fc $\epsilon$ RI $\alpha$  (EqFc $\epsilon$ RI $\alpha$ ). Cultures expressing recombinant EqFc $\epsilon$ RI $\alpha$  were produced as described in US Patent 6,057,127, and the supernatants collected and adjusted to a conductivity of 2.2 mv and a pH of 6.0. The sample was applied to a SP-Sepharose Hi-Trap column (Pharmacia, Peapack, NJ) which had been pre-equilibrated with Buffer A (20 mM Na-PO<sub>4</sub>, pH 6.0). Following application of the sample, the column was washed with Buffer A and the bound proteins eluted using a step-gradient of 0.1M, 0.15M, 0.2M, 0.3M, 0.4M, 0.5M and 1.0M NaCl in Buffer A. Fractions containing EqFc $\epsilon$ RI $\alpha$ , which eluted at 0.2M, were pooled, concentrated, and applied to a Superdex 75 10/30 column which had been equilibrated in PBS. 0.5 ml fractions were collected and analyzed using a Novex 14% SDS-PAGE gel (Invitrogen, Carlsbad, CA) according to manufacturers instructions and fractions containing EqFc $\epsilon$ RI $\alpha$  were pooled.

### Example 8

This example discloses the ability of anti-equine IgE monoclonal antibody H-357 to recognize antigen-bound equine IgE using a solid-phase ELISA. mAb H-357 was produced using the procedure outlined in Example 2 with the exception the mice were immunized with equine IgE instead of canine IgE. H-357 was subsequently shown to be specific for equine IgE. The ELISA's were performed using the following protocol:

The wells of a microtiter plate were coated with the indicated allergen and the plate incubated overnight at 4°C. (All allergens were obtained from Greer Labs, Lenoir, N.C. and coated at 30 PNU/well with the exception of D. farinae which was obtained from Center Labs and coated at 30 allergy units (AU)/well). The following day, excess fluid was removed, 200  $\mu$ l of Assay buffer added to each well and the plate incubated at room temperature (RT) for one hour. The plate was washed four times using Wash buffer, and

100 µl of equine sera (diluted 1:20 in Assay buffer) added to each well and the plate incubated ON at 4°C. The following day, the plate was washed four times using Wash buffer and either biotinylated anti-equine IgE mAb H357 (10 ng), biotinylated, recombinant human Fc epsilon receptor I alpha chain (huFcεRIα) (13 ng) (prepared as described in US Patent 5,945,294) or biotinylated recombinant equine Fc epsilon receptor alpha (EqFcεRIα) (13 ng) (prepared as described in Example 7) was added to the wells. (All molecules were biotinylated following standard procedures using biotin-LC-hydrazide). The plate was incubated at RT for two hours and then washed four times using Wash buffer. Alkaline phosphatase (AP)-conjugated streptavidin (50 ng/well) (KPL Labs) was added to each well and the plate incubated at RT for one hour. Control wells were prepared as described above but without the addition of anti-equine IgE mAb or IgE receptors. The plate was washed using Wash buffer, 100 µl of p-nitrophenyl phosphate (Moss, Inc., Pasadena, MD) added to each well and the plate incubated at RT for one hour. The color reaction was stopped by the addition of 50 µl of 20 mM L-cysteine to each well and the absorbance was read at 405 nM using an ELISA plate reader. The results, with the background values subtracted, are shown in Table 5.

Table 5.

Horse #1			
Allergen	H-357	HuFcεRIα	EqFcεRIα
Black Fly	970	531	575
Horse Fly	1600	695	1134
D. farinae	4035	3807	3942
Horse #2			
Allergen	H-357	HuFcεRIα	EqFcεRIα
Alfalfa	691	448	526
Black Ant	2377	1418	1539
Horse Fly	1484	732	1067
Horse #3			
Allergen	H-357	HuFcεRIα	EqFcεRIα
Black Ant	2138	1719	1972
Horse Fly	450	723	913
Horse #4			
Allergen	H-357	HuFcεRIα	EqFcεRIα
D. farinae	732	414	405
Horse # 5			
Allergen	H-357	HuFcεRIα	EqFcεRIα
Beech	364	388	351
Horse #6			



Allergen	H-357	HuFceR1 $\alpha$	EqFceR1 $\alpha$
Beech	143	140	108
Horse #7			
Allergen	H-357	HuFceR1 $\alpha$	EqFceR1 $\alpha$
Beech	305	345	266

The data demonstrates the absorbance obtained when H-357 is used to detect IgE is similar to that observed when either human or equine FceR1 $\alpha$  is used.

5

### Example 9

This example discloses the specificity of anti-equine IgE monoclonal antibody H-357 for horse IgE. The wells of a microtiter plate were coated, at 50ng/well, with either equine IgE, canine IgE, feline IgE or human IgE in CBC buffer and the plate incubated overnight at 4°C. Excess fluid was removed, 200  $\mu$ l of Assay buffer added to each well and the plate incubated at room temperature (RT) for one hour. The plate was washed four times using Wash buffer and 50  $\mu$ l of serially diluted (from 400 to 3.1 ng/ml in Assay buffer) biotinylated anti-equine IgE mAb (biotinylated by standard procedures using biotin-LC-hydrazide) were added to each well. The plate was incubated at RT for two hours and then washed four times using Wash buffer. HRP-conjugated streptavidin (KPL Labs) was added to each well (10 ng/well) and the plate incubated at RT for one hour. Control wells were prepared as described above but without the addition of anti-equine IgE mAb. The plate was washed using Wash buffer, 100  $\mu$ l of TMB Substrate System added to each well and the plate incubated at RT for 5 minutes. The color reaction was stopped by the addition of 100  $\mu$ l 1M H<sub>3</sub>PO<sub>4</sub> to each well and the absorbance was read at 450 nM using an ELISA plate reader and the results are shown in Table 6. Background levels have been subtracted from the values shown and all values over 100 are considered positive.

Table 6.

mAb H-357 (ng/well)	Coating Immunoglobulin			
	Equine IgE	Canine IgE	Feline IgE	Human IgE
20	2639	208	3	4
10	2567	180	6	0
5	2496	161	0	0
2.5	2077	134	0	0
1.25	1419	98	0	0
.63	848	61	0	0

.31	484	35	3	2
.16	278	15	1	1

This example demonstrates mAb H-357 has a much greater specificity for equine IgE than for canine, feline or human IgE.

#### 5 Example 10

This example discloses the ability of huFceRI $\alpha$  to bind surface-bound IgE from various species. This assay was performed using the basic protocol outlined in Example 9 with the exception that biotinylated-huFceRI $\alpha$ , instead of anti-equine IgE mAb, was used to detect surface-bound IgE. In addition, each well was coated with 100 ng of IgE instead of 50ng IgE, concentrations of receptor were changed and the incubation time with the huFceRI $\alpha$  shortened to one hour. Background levels were measured by preparing control wells with all reagents except huFceRI $\alpha$ . The results of this assay are shown in Table 7. Background values have been subtracted from those shown in the Table.

Table 7.

FceRI $\alpha$ (ng/ml)	Coating Immunoglobulin			
	Human IgE	Canine IgE	Feline IgE	Equine IgE
1300	1482	2286	2973	2123
650	2109	2409	3038	2108
325	1933	2327	3007	2128
162.5	1896	2321	2930	2058
81.3	2128	2307	2741	1867
40.6	1054	2100	2732	1698
20.3	1039	1964	2462	1344
10.2	1462	1911	2443	1132
5.1	1276	1509	2147	881
2.5	783	1050	1447	556
1.3	574	710	891	330

The data demonstrates huFceRI $\alpha$  is capable of recognizing IgE from humans, canines, felines and equines with equal affinity.

#### Example 11

This example discloses the isotypes of the various monoclonal antibodies of the present invention. Isotypes were determined using the IsoStrip Kit (mouse monoclonal

antibody isotyping kit) (Roche Diagnostics Corporation, Indianapolis, IN) according to the manufacturer's instructions. The results of these tests are shown below in Table 8.

Table 8

<b>mAb Designation</b>	<b>Isotype</b>
H99	IgG <sub>1</sub> , k
H100	IgG <sub>1</sub> , k
H101	IgG <sub>1</sub> , k
H102	IgG <sub>1</sub> , k
H103	IgM, k
H104	IgG <sub>1</sub> , k
H106	IgG <sub>1</sub> , k
H107	IgG <sub>1</sub> , k
H189	IgG <sub>2a</sub> , k
H190	IgG <sub>1</sub> , k
H191	IgG <sub>1</sub> , k
H192	IgG <sub>1</sub> , k
H193	IgG <sub>2a</sub> , k
H195	IgG <sub>1</sub> , k
H197	IgG <sub>1</sub> , k
H198	IgG <sub>1</sub> , k
H199	IgG <sub>2b</sub> , k
H200	IgG <sub>1</sub> , k
H201	IgG <sub>1</sub> , k
H202	IgG <sub>1</sub> , k
H203	IgG <sub>1</sub> , k
H204	IgG <sub>1</sub> , $\lambda$
H205	IgG <sub>2b</sub> , k
H206	IgG <sub>2a</sub> , k
H207	IgG <sub>1</sub> , k
H208	IgG <sub>2b</sub> , k
H219 (H187)	IgG <sub>1</sub> , k
H220 (H188)	IgG <sub>2a</sub> , k

Example 12

This example discloses the ability of EqFceRI $\alpha$  to bind surface-bound IgE from various species. This assay was performed using the basic protocol outlined in Example 9 with the exception that biotinylated-EqFceRI $\alpha$ , instead of anti-equine IgE mAb, was used to detect surface-bound IgE. In addition, each well was coated with 100 ng of IgE instead of 50ng IgE, the concentrations of receptor were changed and the incubation time with the EqFceRI $\alpha$  shortened to one hour. Background levels were measured by preparing control

wells with all reagents except EqFceRI $\alpha$ . The results of this assay are shown in Table 8.

Background values have been subtracted from those shown in the Table 9.

Table 9.

FceRI $\alpha$ (ng/ml)	Coating Immunoglobulin			
	Human IgE	Canine IgE	Feline IgE	Equine IgE
1300	142	1805	2612	1838
650	173	1706	2663	1777
325	214	1733	2807	1834
162.5	112	1562	2669	1879
81.3	79	1657	2567	1813
40.6	34	1464	2615	1706
20.3	34	1115	2328	1322
10.2	41	903	2129	1151
5.1	36	599	1653	957
2.5	44	367	1053	684
1.3	7	222	508	385

The data demonstrates EqFceRI $\alpha$  binds canine, feline and equine IgE with high affinity but has a much lower affinity for human IgE.

### Example 13

This example provides the general protocol for performing the immunoassays of the present invention.

Immediately prior to use, the device is wetted with two drops of wash buffer (TRIS 6.05 g/l; NaCl 8.5 g/l; NaN<sub>3</sub> 20 mg/l; Bovine Serum Albumin (BSA) 2.5 g/l; Tween 20 5.0 ml/l; and pH adjusted with concentrated HCl to  $7.5 \pm 0.1$ ). The two drops are added one right after the other, and the drops are allowed to absorb before proceeding.

Two drops of sample (undiluted) are then added, one immediately after the other, and the drops are allowed to absorb before proceeding. Two drops of wash buffer are then added, one at a time, and the drops are allowed to absorb before proceeding.

Two drops of biotinylated mouse anti-canine IgE are then added, one immediately after the other, and the drops are allowed to absorb before proceeding. Two drops of wash buffer are then added, one at a time, and the drops are allowed to absorb before proceeding.

Two drops of streptavidin-alkaline phosphatase are then added, one immediately after the other, and the drops are allowed to absorb before proceeding. Two drops of wash buffer are then added, one at a time, and the drops are allowed to absorb before proceeding.

Two drops of enzyme substrate are then added, one immediately after the other, and the drops are allowed to absorb before proceeding. Wait a minimum of 30 seconds, and add four drops of wash buffer are then added, one immediately after the other. Results can be read immediately and are stable for two hours.

In an alternative embodiment of the invention, the detector reagent is another indicator reagent such as immunogold, or a fluorescent or radioactive compound, and the method of reading the results will change appropriately.

Although the invention has been described with reference to the presently preferred embodiments, it should be understood to those skilled in the art that various modifications can be made without departing from the spirit of the invention. Accordingly, the invention is limited only by the following claims.